

### Technical Data Sheet



Medical Examination
Disposable Nitrile
Powder-Free Gloves

# Nitrile Examination GLOVES

Contact Information:

MANOS MEDICAL SDN. BHD.

Unit 18-3 KL Trillion Office Tower 338 Jalan Tun Razak 50400 Kuala Lumpur Malaysia












## Carton: Gloves Packaging Measurement



#### Product Measurement





#### Product Overview

This product is commonly used in laboratories, hospitals and health centers worldwide. Its balance of protection and tactile sensitivity convinced users at all different kind of applications. While it can be used as a medical examination and category III protective glove, this latex-free protection comes in handy, while you work on glass equipment in laboratories or kitchens. The finger-texturing allows optimal performance by a high a degree of sense of touch throughout the glove film. The finished surface allows grabbing of any smooth and wet surface without losing grip to it, while protecting from any surroundings that should not come in contact with bare skin.

#### Product Specifications

Material	100% Nitrile butadiene rubber (NBR)*		
Sterility	Non-sterile		
Color	Blue		
Package Size	100 pieces / Box		
Weight	Nominal 4.0 grams		
Inner Treatment	Chlorination		
Design & Features	Powder-Free Latex free Ambidextrous Beaded cuff Disposable CE Compliant Figure textured Single use only		
AQL Up to	1.5		
Origin	Malaysia		
Storage	The gloves shall maintain their properties when stored in a dry condition. Avoid direct sunlight		
Shelf-life	5 years from the date of manufacturing		


<sup>\*</sup>This product is not made with Natural Rubber Latex

#### **Quality Standards**

Conforms to ASTM D6319 and EN455 Standards
Conforms to MDD 93 / 42 / EEC Class I and PPE Regulation 2016 / 425 Cat III
Manufactured under QSR (GMP), ISO 9001:2015 and ISO 13485:2016 Quality Management System

#### Glove Sizes

Small, Medium, Large, Extra Large



Size of gloves shall be marked in the check box on the Master Carton with black ink



#### **Physical Dimensions**

Dimensions	Standards			
	GLOVES	ASTM D3578	EN 455	
Length (mm)	Min 230, Min 240, Min 300 ± 10	Min 220 (XS, S) Min 230 (M, L, XL)	Min 240	
Palm Width (mm) Small Medium Large Extra-Large	84 ± 3 94 ± 3 105 ± 3 113 ± 3	80 ± 10 95 ± 10 110 ± 10 120 ± 10	80 ± 10 95 ± 10 110 ± 10 ≥ 110	
Thickness: Single Wall (mm) Finger Palm	Min 0.05 Min 0.05	Min 0.05 Min 0.05	N/A N/A	

#### **Physical Properties**

Property	GLOVES	ASTM D6319	EN 455
Elongation at Break (%) Before Aging After Aging	Min 500	Min 500	N/A
	Min 400	Min 400	N/A
Tensile Strength (MPa) Before Aging After Aging	Min 14	Min 14	N/A
	Min 14	Min 14	N/A
Median Force at Break (N) Before Aging After Aging	N/A	N/A	Min 6
	N/A	N/A	Min 6

#### Instructions

Non-Sterile; Not Surgical Gloves

When wearing gloves, do not wear rings or other accessories, pay attention to trim nails Wash hands before and after use

NOTE: Not for use in diagnostic procedures.

#### Symbols



Disclaimer: JKK GREENGOLD SDN. BHD. believes that the information in this Technical Data Sheet is accurate and represents the best and most current information available to us. JKK GREENGOLD SDN. BHD. makes no representations or warranties either express or implied, regarding the suitability of the material for any purpose or the accuracy of the information contained within this document. Accordingly, JKK GREENGOLD SDN. BHD. shall not be responsible for damages resulting from use of or reliance upon this information.